# \*Strep B Carrot Broth™

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1. Regan, J.A., Klebanoff, M.A., Nugent, R.P. 1991. *The epidemiology of group B streptococcal colonization in pregnancy*. Vaginal Infections and Pregnancy Study Group. *Obstet. Gynecol.*; 77:604-10.
2. Schrag, S.J., E.R. Zell, R. Lynfield, et al. 2002. *A population-based comparison of strategies to prevent early-onset group B streptococcal disease in neonates*. *N. Engl. J. Med.* 25;347(4):233-9.
3. Schuchat, A. 2001. *Group B streptococcal disease: from trials and tribulations to triumph and trepidation*. *Clin. Infect. Dis.* 5;33(6):751-6.
4. The Centers for Disease Control and Prevention. 2010. *Prevention of Perinatal Group B Streptococcal Disease*. Revised Guidelines. MMWR 59 (RR-10). Internet: <http://www.cdc.gov/mmwr/pdf/rr/rr5910.pdf>
5. Overman, S.B., D.D. Eley, B.E. Jacobs, J.A. Ribes. 2002. *Evaluation of methods to increase the sensitivity and timeliness of detection of Streptococcus agalactiae in pregnant women*. *J. Clin. Microbiol.*; 40(11):4329-31.
6. de la Rosa, M., M. Perez, C. Carazo, L. Pareja, J.I. Peis and F. Hernandez. 1992. *New Granada Medium for detection and identification of group B streptococci*. *J. Clin. Microbiol.*; 30:1019-1021.
7. Garcia Gil, E., M.C. Rodriguez, R. Bartolome, B. Berjano, L. Cabero and A. Andreu. 1999. *Evaluation of the Granada Agar plate for detection of vaginal and rectal group B streptococci in pregnant women*. *J. Clin. Microbiol.*; 37:2648-2651.
8. Rosa-Fraile, Manuel, J. Rodriguez-Granger, M. Cueto-Lopez, A. Sampedro, E. Biel Gaye, J.M. Haro and A. Andreu. 1999. *Use of Granada Medium to detect group B streptococcal colonization in pregnant women*. *J. Clin. Microbiol.*; 37:2674-2677.
9. Rosa-Fraile, Manuel, A. Sampedro, J. Varela, M. Garcia-Pena, and G. Gimenez-Gallego. 1999. *Identification of a peptide pigment from mammal albumins responsible for enhanced pigment production by group B streptococci*. *Clin. Diag. Lab. Imm.*; 6:425-426.
10. Islam, AKMS. 1977. *Rapid recognition of group B streptococci*. *Lancet* 309(8005):256-257.
11. B. Spellerberg, B. Pohl, G. Haase, S. Martin, J. Weber-Heynemann and R. Luticken. 1999. *Identification of Genetic Determinants for the Hemolytic Activity of Streptococcus agalactiae by ISS1 Transposition*. *J. Bacteriol.*; 181: 3212-3219.
12. Hardy Diagnostics and Fleury Medical Diagnostic Center. 2004. *Evaluation of Three Methods for Recovery of Group B Streptococci*. A poster presentation at American Society for Microbiology 104th General Meeting, New Orleans, LA.
13. National Center for Infectious Diseases, et al. 2002. Revised Guidelines from CDC: *Prevention of Perinatal Group B Streptococcal Disease*. MMWR Recommendations and Reports/51(RR11);1-22. Internet: <http://www.cdc.gov/mmwr/preview/mmwrhtml/rr5111a1.htm>.
14. Rosa-Fraile, M. et al. 2005. *Specimen Storage in Transport Medium and Detection of Group B Streptococci by Culture*. *J. Clin. Microbiol.*; 43:928-930.
15. Schreckenberger et al. 2005. *Evaluation of Strep B Carrot Broth™ and LIM Broth Methods for Recovery of Group B Streptococci (GBS): Results of a Multi-Center Trial*. A poster presentation at American Society for Microbiology, 105th General Meeting, Atlanta, GA.
16. DiPersio, J. 2005. *GBSPreservation by Strep B Carrot Broth™ Inoculated with Patient Specimens*. Unpublished data.
17. Facklam, R. et al. 2006. *Evaluation of Accuracy of Strep B Carrot Broth™ in the Detection of Different Serotypes of Group B Streptococci (GBS)*. A poster presentation at American Society for Microbiology, 106th General Meeting, Orlando, FL.
18. Young, Uh et al. 1998. *Serotypes and Biochemical Reaction Patterns of Group B Streptococci*. *Korean J. Clin. Path.*; 18:386-390.
19. Merrit, K. and Jacobs, N. 1978. *Characterization and Incidence of Pigment Production by Human Clinical Group B Streptococci*. *J. Clin. Microbiol.*; Vol. 8, No. 1, p. 105-107.
20. Noble, M., Bent, J., West, A. 1983. *Detection and identification of group B streptococci by use of pigment production*. *J. Clin. Path.*; 36:350-352.
21. de la Rosa, M., R. Villareal, D. Vega, C. Miranda, and A. Martinezbrocal. 1983. *Granada Medium for Detection and Identification of Group B Streptococci*. *J. Clin. Microbiol.*; Vol. 18, No. 4, p.779-785.
22. Rosa-Fraile, Manuel, J. Rodriguez-Granger, A. Haidour-Benamin, J.M. Cuerva, and A. Sampedro. 2006. *Granadaene: Proposed Structure of the Group B Streptococcus Polyenic Pigment*. *J. Clin. Microbiol.*; Vol. 72, No. 9, p.6367-6370.
23. da Gloria Carvalho, M., R. Facklam, D. Jackson, B. Beall, and L. McGee. 2009. *Evaluation of Three Commercial Broth Media for Pigment Detection and Identification of a Group B Streptococcus (Streptococcus agalactiae)*. *J. Clin. Microbiol.*; Vol. 47, No. 12, p.4161-4163.
24. Block, T., E. Munson, A. Culver, K. Vaughan, and J. Hryciuk. 2008. *Comparison of Carrot Broth- and Selective Todd-Hewitt Broth-Enhanced PCR Protocols for Real-Time Detection of Streptococcus agalactiae in Prenatal Vaginal/Anorectal Specimens*. *J. Clin. Microbiol.*; Vol. 46, No. 11, p.3615-3620.
25. Church, D.L., H. Baxter, T. Lloyd, B. Miller, and S. Elsayed. 2008. *Evaluation of Strep B Carrot Broth™ versus Lim Broth for Detection of Group B Streptococcus Colonization Status of Near-Term Pregnant Women*. *J. Clin. Microbiol.*; Vol. 46, No. 8, p.2780-2782.
26. Czerepuszko, D.J., and M.J. Lewis. 2010. *Comparison of LIMBroth with PNA FISH to Carrot Broth with PNA FISH for Identification of Group B Streptococcus in Prenatal Vaginal/Rectal Specimens*. A poster presentation at American Society for Microbiology, San Diego, CA.

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Copyright © Hardy Diagnostics. All rights reserved.

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Copyright © Shimadzu Diagnostics Corporation. All rights reserved.

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